

**Assignment 4 – Due Friday, 3/09/12, noon.****Thermodynamic Analysis of RNaseA Denaturation by UV-Vis Difference Absorption Spectroscopy (and Differential Scanning Calorimetry).**

The accompanying excel file (Origin\_Assign\_4\_Data.xlsx) contains two sets of data. The first set of data represents the change in molar absorptivity at 280 nm ( $\Delta\epsilon_{280}$ ) versus temperature of ribonuclease A (RNaseA). The second data set are the results of a DSC experiment in which the heat capacity ( $C_p$ ) was measured as a function of temperature.

**1. Use the Data for RNase to estimate:**

- a.  $T_M$ , the transition temperature (where  $K = [D]/[N] = 1$ ),
- b.  $\Delta H$  and  $\Delta S$  of denaturation (at the transition temperature), and
- c.  $\Delta C_p$  of denaturation.
- d. Provide a write-up that outlines the procedure you used to estimate these parameters. Include (d1) your rationale and procedure for calculating baselines, (d2) the method for determining  $K$  as a function of temperature, (d3) a plot of  $\Delta G$  as a function of temperature, and (d4) the procedure for determining  $\Delta H$  and  $\Delta S$  of denaturation (at the transition temperature).
- e. Calculate the temperature of maximum stability.
- f. Generate a protein stability plot ( $\Delta G$  versus  $T$ ) based on the fit parameters generated in part d and use the plot to predict the  $T_M$  of cold denaturation.

(Assume two-state behavior for RNaseA)

**Extra Credit.** Carry out an analysis of DSC data for RNase to determine,  $T_M$ ,  $\Delta H_M$ ,  $\Delta S_M$ ,  $\Delta C_p$ , the temperature of maximum stability, and the  $T_M$  of cold denaturation ( $T'_G$ ) in that manner that you did for problem 1. In addition to the vant Hoff enthalpy change, calculate  $\Delta H_{cal}$ , the calorimetric change in enthalpy.

**Conditions of the DSC Experiment**

Scan Rate = 1 deg per minute

Cell Volume = 1.411 mL

[RNase] = 63  $\mu$ M